Full Blade Master Blade Total Blade Cyl Burr Ball Burr

Puntas de Shaver

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Instructions for use
ENGLISH

M-350-00845 (Rev. 03) / 20.04.2024

Fergus

ESPAÑOL

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ENGLISH

1. DEVICE DESCRIPTION

Fergus shaver blades and burrs are specifically designed with different shaver configurations and cutting windows for soft tissue, cartilage and bone resection during arthroscopic surgical procedures.

The Fergus shaver blades and burrs include the Full Blade, Master Blade, Total Blade, Cyl Burr and Ball Burr models.

The **Full Blade** model is designed for the meniscal, cartilage, and synovial membrane and plicae debridement.

The **Master Blade** model is designed for joint soft tissue debridement.

The **Total Blade** model is designed for articular cartilage debridement as well as for the resection of the periosteum of the acromion.

The **Cyl Burr** model is designed for swift and aggressive resection of the femoral bone in intercondylar notch remodeling, or of the distal clavicle, and for subacromial decompression procedures.

The **Ball Burr** model is designed to be used for bone resection of the intercondylar notch and subacromial decompression.

2. PRODUCT PRESENTATION

The devices are packaged in blisters which are sterilized using ethylene oxide.

Check the model and the diameter of the shaver blade/burr in each of the individual device labels.

3. INDICATIONS FOR USE

The Fergus shaver blades and burrs are designed for soft and bone tissue resection during arthroscopic surgical procedures.

4. CONTRAINDICATIONS

Insufficient bone quality or quantity: the physician shall carefully assess the bone quality of the patients with bone immaturity before starting the orthopedic surgical procedure. The growth cartilage must not be deviated, interrupted or altered.

Limited blood irrigation and/or previous infection which may delay consolidation.

5. WARNINGS

Read the instructions for use carefully and completely before using the product.

The device must be used only by a trained medical professional.

The use of this device is restricted to physicians or people authorized by a physician.

DO NOT use the device if the packaging is open or damaged since product sterility will be lost.

Shaver blades and burrs have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their efficiency and increase the risk of inappropriate sterilization and cross-contamination.

The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative or governmental requirements.

Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.

The product must be stored in the original packaging.

Bio-hazardous waste and contaminated surgical materials must be disposed of safely and as per the policies of the site where the surgical procedure is performed.

The manufacturer shall not be liable for any effects or consequences which may arise from deviating from these specific instructions.

6. PRECAUTIONS

Surgeons must follow their professional judgement to determine the appropriate device model to use for the surgery.

Device components must not be interchanged.

Do not allow any rotating parts of the device to come into contact with metal parts, since this could cause irreversible damage to both objects. Contact between metallic parts could change the sharpness of the cutting edge of the blade and/or cause the tip to break inside the body.

Do not apply excessive force to the device.

Do not use the device without irrigation flow. Otherwise, you will increase the temperature of the device and, therefore, cause burns to the patient.

Do not handle the cutting terminal end.

7. SUPPLY AND STORAGE

The devices are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide. In order to maintain product quality and effectiveness it is supplied in a blister. After the sterilization, the product is protected during transportation and storage by placing it inside a cardboard box. Labels may be added for identification.

Keep the products away from any sources of heat or humidity.

Room temperature storage is advised for this product.

IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.

DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.

Product shelf life is 3 years.

8. POSSIBLE COMPLICATIONS

Additional damage due to applying excessive force while using the device.

Infection, both deep and superficial.

9. MODE OF USE

Open the packaging, remove the sterile product carefully, and place it in a sterile area to check its integrity.

Mount the shaver blade/burr to the handpiece following these instructions:

a. Coupling:

To place the shaver blade/burr in the handpiece, it must be inserted so that the locking structure is aligned with the handpiece cavity. If it does not become locked in place successfully, remove the shaver blade/burr, turn its internal shaft to align it with the engine output shaft and place it again, making sure it locks securely.

NOTE: For optimal performance during surgery, we recommend using shaver blades at 3000 rpm and shaver burrs at 6000 rpm.

b. Uncoupling:

To uncouple the device from the handpiece, press the button on the handpiece.


10. MATERIAL SPECIFICATIONS

Fergus shaver blades and burrs are made with stainless steel, polycarbonate, aluminum and polytetrafluoroethylene.

11. LEARN MORE

To learn more about this product, visit our website: www.fergusmedical.com

12. SYMBOLS USED IN LABELS

	CATALOG NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	PRECAUTION
	EXPIRATION DATE
	DO NOT REUSE
	MANUFACTURING DATE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE. SINGLE STERILE BARRIER SYSTEM
	MANUFACTURER
	CHECK INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	DO NOT RE-STERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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PROMEDON S.A.