# **Knotix** PK **QuickFit** PK **Knotix** SPD PK **SpinFit** PK

---

Anclas de PEEK

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Instructions for use  
ENGLISH

M-351-00148(Rev.03) / 10.01.2023

Fergus

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## 1. DESCRIPTION OF THE DEVICE

The Fergus-brand PEEK anchors include the Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK models. These are devices that are designed with threads or tabs (i. e., non-threaded) and with eyelets to allow for the passage of suture. They are supplied sterile and inside a disposable insertion device to aid in the placement of the implant.

Some PEEK anchor models may include non-absorbable suture, while others may not include sutures in the kit and feature a suture threader to help threading a suture afterwards.

The PEEK anchors are designed to attach soft tissue to the bone.

The **Knotix PK** model features a threaded anchor and an internal eyelet for passing two CordFit #2 sutures (USP standard). Sutures may come in blue, white or white-blue colors. The four tips of these two sutures are connected to the suture retention device located in the handle of the insertion device. This model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 5.5 mm diameter, 14.7 mm length.

Figure 1: Knotix PK PEEK anchor.

The **Knotix SPD PK** model comes in two variants:

- a non-threaded anchor with tabs and a non-absorbable CordFit #2 (USP standard) suture threaded through the external eyelet of the anchor. Sutures may come in blue, white or white-blue colors. The 2 tips of the suture are connected to the suture retention device located in the handle of the insertion device. The Knotix SPD PK model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 3.0 mm diameter, 12.4 mm length.
- a non-threaded anchor with tabs and an internal eyelet for passing two non-absorbable CordFit #2 sutures (USP standard). Sutures may come in blue, white or white-blue colors. The four tips of these two sutures are connected to the suture retention device located in the handle of the insertion device. The Knotix SPD PK model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 3.0 mm diameter, 12.2 mm length.

Figure 2: Knotix SPD PK PEEK anchor  
(top: one suture, bottom: two sutures)

The **QuickFit PK** PEEK anchor is also a two-part implant: a non-threaded tabbed anchor and the eyelet to pass the non-absorbable suture. The suture is not included in the kit and must be threaded by the physician. This model includes a suture threader built into the eyelet to aid in threading the suture. The eyelet is attached to the insertion device by means of a small thread in the rod. During implantation, the anchor is connected to the eyelet, which allows for the suture to remain in place due to an interference phenomenon between the bone tissue and the implant, with no knots required. This model is available in the following dimensions: 4.5 mm diameter and 24 mm length, or 2.9 mm diameter and 15.5 mm length. Dimensions are identified in product labels.

Figure 3: QuickFit PK PEEK anchor

The **SpinFit PK** PEEK anchor is a two-part implant: a threaded anchor and an eyelet to pass the suture. The suture may be included in the kit or not. If not, the physician must thread it in the eyelet. This model is available in the following dimensions: 4.7 mm diameter and 19.1 mm length, or 3.5 mm diameter and 15.8 mm length. Dimensions are identified in product labels.

The anchor eyelet is distal, separated from the anchor and attached to the insertion device by means of an eyelet retention suture, for the 4.7 mm diameter anchor, or a thread, for the 3.5 mm diameter anchor.

The SpinFit PK anchor model that does not include a suture, comes with a suture threader to aid the physician in passing the sutures through the eyelet. If the model includes a suture, it is a TapeFit Max, which is a combination of a thread-type suture on the ends and a tape-type suture in the center. The width in the center of the suture is 2 mm, and the length of the suture is 30". It is available in white-blue and white-black colors.

During implantation, the anchor moves and is attached to the eyelet when they come into contact. This causes the suture to become fixed in place by means of an interference phenomenon between the bone tissue and the implant thread, with no knots required.

Figure 4: SpinFit PK PEEK anchor with suture

Figure 5: SpinFit PK PEEK anchor without suture

## 2. PRODUCT PRESENTATION

Contents:

- **Knotix PK models:** one implantable PEEK anchor with two CordFit #2 implantable sutures and one insertion device.
- **Knotix SPD PK models:**
  - One implantable PEEK anchor with one CordFit #2 implantable suture and an insertion device.
  - One implantable PEEK anchor with two CordFit #2 implantable sutures and one insertion device.
- **QuickFit PK models:** one implantable PEEK anchor, one insertion device and one suture threader.
- **SpinFit PK model:**
  - One implantable PEEK anchor, one insertion device and one suture threader.
  - One implantable PEEK anchor with one TapeFit Max suture and one insertion device.

Check the individual labels of the device to find the model, anchor size, materials, whether the suture or the insertion device are included or not, and the color of the suture. All components are provided sterile and ready for use.

## 3. INDICATIONS FOR USE

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK PEEK anchors are intended for attaching soft tissue to the bone. They are designed for the following surgical procedures:

- **Shoulder joint:**
  - Rotator cuff repair,

- Glenoid labrum SLAP lesion repair, and
  - Glenoid labrum Bankart lesion repair.
- **Ankle joint:**
- Lateral instability repair,
  - Medial instability repair,
  - Instability repair augmentation, and
  - Achilles tendon repair.

Surgeons must follow their professional judgment to determine the appropriate implant size and model based on the specific indication and the history of the patient. However, the 4.5 mm diameter Knotix PK and QuickFit PK or the 4.7 mm diameter SpinFit PK PEEK anchor models are recommended for rotator cuff repair (since, due to their size, their properties are a better fit for attaching these tendons).

For SLAP and Bankart lesion repairs, the 2.9 mm diameter Knotix SPD PK and QuickFit PK, or the 3.5 mm SpinFit PK, PEEK anchors are advised (due to their smaller size, their implementation in the glenoid cavity is easier).

For lateral and medial ankle instability repairs, the Knotix PK and Knotix SPD PK are recommended. For Achilles tendon repair, all anchor models can be used, except for the 2.9 mm diameter QuickFit PK models and the 3 mm diameter Knotix SPD PK models. Finally, for ankle instability repair augmentations, the SpinFit PK and the QuickFit PK anchor models, in both diameters, are advised.

The anchor creates an anchoring point for the suture in the bone surface, while the sutures attach the intervening tendons to the anchoring point.

The insertion device allows for the implantation of the anchor and the suture in the right place.

#### 4. CONTRAINDICATIONS

The device should not be implanted in the following cases:

- Insufficient bone quality or quantity for attachment. The physician must evaluate carefully the bone quality of patients with bone immaturity before starting with the orthopedic surgical procedure, and must avoid deviating, interrupting or altering the growth cartilage.
- Tissue quality or pathologies of the soft tissue to be reconstructed which may affect the safety of the attachment by means of suture.
- Limited blood irrigation or previous infection which may delay healing.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions of their physical activity.
- Surgical procedures not mentioned in the "Indications for Use" section.
- Inappropriate neuromuscular conditions.

## 5. WARNINGS

- Read the instructions for use carefully and completely before use.
- The device must be used only by a medical professional trained in the surgical techniques indicated for placing these implants.
- The law restricts the sale of this device to physicians or people authorized by a physician.
- DO NOT use the product if the packaging is open or damaged since product sterility may be lost.
- The products have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative and/or governmental requirements.
- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- Make sure that the product is disposed as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- The product must be stored in the original packaging.
- Patient sensitivity to device materials must be considered before implantation.
- All surgical instruments used in the procedure to install this implant must be sterile and suitable for use in aseptic environments.
- Product presentation must not be altered in a way different from what is stated in the surgical technique, since this could affect implant performance.
- If the implant detaches from the insertion device, DO NOT REATTACH, since this could impair the effectiveness and increase the risk of sterility loss.
- The fixation provided by the device must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing the device and the bone.
- The patient will receive detailed instructions about the use and limitations of the device.
- If the device is to be explanted, the risks of the possibility of a second procedure on the patient must be considered. Device extraction must be accompanied by an appropriate surgical and post-surgical approach.
- Opened devices and implants must be disposed of, even if they have not been used or handled, since this could increase the risk of inappropriate sterilization and cross-contamination.
- Bio-hazardous waste, such as explanted devices, sutures, needles and contaminated surgical materials, must be disposed of safely and as per site policy.


- The anchor and/or the insertion device may break when the insertion site is not properly prepared with the appropriate instrument before using the implant.
- Do not use the product if any faults are detected.
- Do not use products or sutures that are not recommended by the manufacturer, since implant performance for the repair cannot be guaranteed in this case.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the manufacturer.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, which may affect patient safety and implant performance.

## 6. PRECAUTIONS


- Surgeons must follow their professional judgment to determine the appropriate implant model, quantity and size based on the specific indication and the history of the patient.
- Surgeons are advised to review the product-specific surgical technique before the procedure.
- PEEK anchor implantation requires preparation in the insertion site.
- The surgeon must consider the anchor size to select the punch, drill bit or bone drill, since an incomplete insertion of the anchor may cause a deficient performance of the anchor or an inappropriate implantation.
- The surgeon must consider the manufacturer-recommended placement/implantation systems in order to implant the device properly.
- If the bone drill or punch is not inserted to the appropriate depth, the suture anchor may become damaged or broken.
- In order to prevent damage on the anchor and/or insertion device, the surgeon must introduce them in the same orientation as the hole drilled in the bone.
- Using excessive force during insertion may cause faults in the suture anchor or the insertion device.
- The surgeon must check that the anchor, the eyelet and the sutures are properly coupled with the insertion device. If not, they may become damaged or broken.
- Do not uncouple the suture until the anchor has been introduced completely into place, since this could affect the implantation.
- For SpinFit PK anchors only:
  - Insert the capillary of the insertion device into the bone cavity up to the laser marking of the capillary or rod. Check and adjust suture tension at this

point, since this value could vary while the anchor advances to its final position and readjustment will not be possible.

- Check that the anchor comes into contact with the bone before introducing the body of the implant in the prepared bone cavity.
- For QuickFit PK anchors only:
  - Before impacting the anchor, remember to remove the safety from the proximal end of the handle.
  - Make sure that, after the implantation impacts, the external capillary has descended fully into the distal transversal laser marking.
  - Insert the internal rod of the insertion device into the bone cavity until the body of the anchor comes into contact with the bone or up to the laser marking of the rod. Check and adjust suture tension at this point, since this value could vary while the anchor advances to its final position and readjustment will not be possible.

- Check that the anchor body comes into contact with the bone before implanting it fully in the prepared bone cavity.
- Sutures may break if handled or controlled using sharp or cutting instruments.
- The decision of removing the implant must be made considering the potential risk that involves a second surgical procedure for the patient. After the implant is removed, appropriate post-operative controls must be carried out.
- Before using the device, examine it for damage. Do not use if damaged.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper device selection and installation, are essential for satisfactory use of the implant.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.

## 7. SUPPLY AND STORAGE

- PEEK anchors are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- In order to maintain product quality and effectiveness, they are marketed in a double Tyvek pouch. After the sterilization, the product is protected from light and heat by placing it inside a cardboard box (external packaging). The product is identified with labels. Additionally, to prevent damage to the Tyvek pouches, a plastic protector is placed to contain the insertion device and cover the implant.
- Keep the product away from any sources of heat or humidity.
- Room temperature storage is advised.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- The implant has a 3 year shelf life.

## 8. MAGNETIC RESONANCE (MR) ENVIRONMENTS

### • MR Safety and Compatibility

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK anchors are manufactured with polyether ether ketone (PEEK), which makes them MR safe.

Sutures are made with UHMWPE or a mixture of UHMWPE with polypropylene or polyamide (nylon), which makes them MR safe.

## 9. POTENTIAL COMPLICATIONS

All surgical procedures involve a certain risk. Some of the complications that could arise are:

- Deep and superficial infections
- Atypical reactions to foreign bodies
- Pain on the site of the incision or the procedure
- Implant loosening and/or detachment
- Bone lesions
- Incomplete or inappropriate lesion correction, or lesion recurrence.

## 10. INFORMATION FOR PATIENTS

The surgeon is responsible for informing the patient and his/her representatives before surgery about the possible complications related to the implantation of the device (see section "Potential complications").

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated

with the implantation and use of the device.

Post-operative care is important. The patient will receive instructions from the physician about the limitations required by the implant, and will be warned about the necessary precautions they must apply for healing.

The patient must be warned that, if explantation of the device is deemed necessary, the risks of the possibility of a second surgical procedure on the patient must be considered. Device extraction must be accompanied by an appropriate post-operative approach.

Performing additional surgeries does not guarantee a complete repair of the complications.

## 11. SURGICAL TECHNIQUE

The surgeon must select the appropriate approach depending on the repair to be performed. For foot and ankle procedures, the physician may opt for an open or arthroscopic approach. For shoulder surgeries, an arthroscopic approach is recommended.

Additionally, the surgeon must select the appropriate number of implants and fixation configuration based on the type of lesion, the medical history of the patient and other criteria that the professional may deem important for a satisfactory use of the product during the surgical procedure.

The following surgical technique instructions refer only to device implantation, not considering other procedures that are critical to perform the surgical procedure, such as patient preparation, working portals, procedure finalization, etc.

### ***Knotix PK***

The Knotix PK anchor is recommended to perform rotator cuff repairs in the shoulder joint, medial and lateral ankle instability repairs, and Achilles tendon repair.

Steps for proper Knotix PK anchor implantation:

- 1)** Identify the anchor insertion site and place the 5.5 mm Knotix PK punch. Impact the end of the handle of the 5.5 mm Knotix PK punch successively until the first laser-marked line of the capillary of the instrument has been inserted completely in order to create the hole in the bone.
- 2)** Position the anchor in the implantation site.
- 3)** Turn the handle of the device clockwise in order to implant the anchor in the bone tissue by means of threading. Applying a slight forward (proximal-distal) pressure is advised while the handle is turned in order to aid in the implantation.

**ATTENTION:** Ensuring a complete anchor insertion is important to prevent possible damage in the surrounding tissue or implant detachment. To achieve this, it must be threaded until the

transversal laser marking located in the distal end of the capillary (of the insertion device) has been completely inserted in the hole where the anchor is held.

**4)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor.

**5)** For arthroscopic approach: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached; tie the knot, and then cut the excess suture (if appropriate).

For open approach: the suture must be threaded in a needle in order to be passed through the tissue, and attached to the anchoring point. Finish the attachment by tying surgical knots.

Note: the needle is not included in the kit.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### ***Knotix SPD PK***

The Knotix SPD PK anchor is intended to perform glenoid labrum repairs for shoulder instability, and medial and lateral ankle instability.

Steps for proper Knotix SPD PK anchor implantation (with one or two sutures):

**1)** Identify the implantation site for the anchor and position the 3.2 mm implantation guide in the selected insertion site.

**2)** Insert the 3.0 mm drill bit for Knotix SPD through the inside of the guide and drill a hole in the insertion point.

**3)** Insert the anchor through the guide and impact the handle (of the insertion device) successively in order to be able to implant the anchor in the bone tissue by means of impact.

ATTENTION: Ensuring a complete anchor insertion is important to prevent possible damage in the surrounding tissue or implant detachment. To achieve this, it must be impacted until the transversal laser marking located in the distal end of the capillary (of the insertion device) has been completely inserted in the hole where the anchor is held.

**4)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor and the implantation guide.

**5)** For arthroscopic approach: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached; tie the knot, and then cut the excess suture.

For open approach: the suture must be threaded in a needle in order to be passed through the tissue, and attached to the anchoring point. Finish the attachment by tying surgical knots.

NOTE: the needle is not included in the kit.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **SpinFit PK** (4.7 mm in diameter)

The 4.7 mm diameter SpinFit PK anchor is recommended to perform rotator cuff repairs in the shoulder joint, ankle instability repair augmentation, and Achilles tendon repair.

Steps for proper SpinFit PK anchor implantation:

**1)** Identify the anchor insertion site and place the 4.7 mm SpinFit PK punch. Impact the end of the handle of the 4.7 mm SpinFit PK punch successively until the second laser-marked line of the capillary of the instrument has been inserted completely in order to create the hole in the bone.

**2)** Remove the 4.7 mm SpinFit PK punch and position the SpinFit PK anchor inside the hole created in the bone until the laser line marked in the anchor capillary is not arthroscopically visible. For SpinFit PK without TapeFit Max suture: Before this step, the suture holding the tissue must be passed through the anchor eyelet using the disposable threader connected to the capillary of the insertion device.

ATTENTION: Check for the desired suture tension, since, from the insertion stage on, the suture will be secured and you will be unable to readjust attachment tension.

**3)** Turn the proximal end of the handle clockwise holding the insertion device to be able to move forward towards the inside of the hole in the bone. Turn until the implantable anchor is not arthroscopically visible and is inserted completely in the hole in the bone.

**4)** Decouple the eyelet retention suture from the proximal end of the handle and pull from one of the ends of the suture to remove it.

**5)** Release the TapeFit Max suture from the handle tabs and remove the insertion device from the anchor.

For SpinFit PK without TapeFit Max suture: cut the excess suture.

**6)** For repairs: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached.

For augmentation: skip the aforementioned step.

**7)** Cut the excess suture.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **SpinFit PK** (3.5 mm in diameter)

The 3.5 mm diameter SpinFit PK anchor is

recommended to perform glenoid labrum repairs in the shoulder, ankle instability repair augmentation, and Achilles tendon repair.

Steps for proper SpinFit PK anchor implantation:

**1)** Identify the insertion site for the anchor and position an anchor insertion guide with an inner diameter of, at least, 3.2 mm, in the insertion site.

**2)** Insert a 2.9 mm drill bit through the inside of the guide and drill a hole 15.8 mm in depth in the insertion point. Remove the guide and the drill bit.

**3)** Thread the suture holding the tissue through the eyelet of the SpinFit PK anchor and hold the sutures in the suture retention area of the handle.

**ATTENTION:** Check for the desired suture tension, since, from the insertion stage on, the suture will be secured and you will be unable to readjust attachment tension.

**4)** Turn the proximal end of the handle clockwise holding the insertion device to be able to move forward towards the inside of the hole in the bone. Turn until the implantable anchor is not arthroscopically visible and is inserted completely in the hole in the bone.

**5)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor by rotating it counter-clockwise until it detaches from the eyelet. Cut the excess suture.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **QuickFit PK** (2.9 mm in diameter)

The 2.9 mm diameter QuickFit PK anchor is intended to perform glenoid labrum repairs for shoulder instability, and ankle instability repair augmentation.

Steps for proper QuickFit PK anchor implantation:

**1)** Identify the insertion site for the anchor and position an anchor insertion guide with an inner diameter of, at least, 3.2 mm, in the insertion site.

**2)** Insert a 2.9 mm drill bit through the inside of the guide and drill a hole 15.5 mm in depth in the insertion point. Remove the guide and the drill bit.

**ATTENTION:** To achieve a 15.5 mm deep perforation, the drill bit must advance until the drill comes into contact with the guide handle.

**3)** Thread the suture holding the tissue through the eyelet of the QuickFit PK anchor and hold them in the suture retention area of the handle.

**ATTENTION:** Check for the desired suture tension, since, from the insertion stage on, the suture will be secured and you will be unable to readjust attachment tension.

**4)** Position the anchor over the perforation made

and impact the handle (of the insertion device) successively to move the end of the inner rod (where the eyelet is located) forward into the perforation. Impacts should continue until the first transversal laser marking (located in the rod) has been inserted completely in the perforation.

**5)** Remove and discard the safety clip located in the proximal end of the handle.

**6)** Impact the handle (of the insertion device) successively to advance it forward and perform the final implantation of the anchor. Impacts should continue until the second transversal laser marking (located in the external capillary) has been inserted completely in the perforation.

**ATTENTION:** In order to prevent possible damage in the surrounding tissue or implant detachment, it is important that you ensure a complete insertion of the anchor. To achieve this, the impacts must be performed ensuring that the transversal laser markings are not visible (as appropriate) after each of the steps above.

**7)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor by rotating it counter-clockwise until it detaches from the eyelet. Cut the excess suture.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **QuickFit PK** (4.5 mm in diameter)

The 4.5 mm diameter QuickFit PK anchor is recommended to perform rotator cuff repairs in the shoulder joint, ankle instability repair augmentation, and Achilles tendon repair.

Steps for proper QuickFit PK anchor implantation:

**1)** Identify the anchor insertion site and place the 4.5 mm QuickFit PK punch. Impact the end of the handle of the 4.5 mm QuickFit PK punch successively until the second laser-marked line of the capillary of the instrument has been inserted completely in order to create the hole in the bone.

**2)** Thread the suture holding the tissue through the eyelet of the QuickFit PK anchor and hold them in the suture retention area of the handle.

**ATTENTION:** Check for the desired suture tension, since, from the insertion stage on, the suture will be secured and you will be unable to readjust attachment tension.

**3)** Position the anchor over the perforation made and impact the handle (of the insertion device) successively in a distal-proximal direction to move the end of the inner rod (where the eyelet is located) forward into the perforation. Impacts should continue until the first transversal laser marking (located in the rod) has been inserted

completely in the perforation.

**4)** Remove and discard the safety clip located in the proximal end of the handle.

**5)** Impact the handle (of the insertion device) successively in a distal-proximal direction to advance it forward and perform the final implantation of the anchor. Impacts should continue until the second transversal laser marking (located in the external capillary) has been inserted completely in the perforation.

**ATTENTION:** In order to prevent possible damage in the surrounding tissue or implant detachment, it is important that you ensure a complete insertion of the anchor. To achieve this, the impacts must be performed ensuring that the transversal laser markings are not visible arthroscopically (as appropriate) after each of the steps above.

**6)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor by rotating it counter-clockwise until it detaches from the eyelet. Cut the excess suture.

After use, discard the product, sutures, needles, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

## 12. MATERIAL SPECIFICATIONS

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK PEEK anchors are manufactured with polyether ether ketone (PEEK), as per the ASTM F2026 standard. CordFit #2 (USP standard) sutures are manufactured with ultra-high-molecular-weight polyethylene (UHMWPE). Sutures used may be blue, white or white-blue. The tint used for blue sutures and the blue part of the white-blue suture is 73.1015 chromium-cobalt-aluminum oxide. White sutures are not tinted.

The 30" TapeFit Max sutures are manufactured with a mixture of ultra-high-molecular-weight polyethylene (UHMWPE) and polypropylene, for white-blue suture, or a mixture of ultra-high-molecular-weight polyethylene (UHMWPE) and polyamide (nylon), for white-black suture. The blue tint used for these sutures is phthalocyanine, while the black tint used is hematein. White fibers in the suture are not tinted.

The insertion device is made up of the following materials:

- Handle: it is manufactured with medical grade ABS, Lustran 348. There are different colors to identify the anchor model for which the insertion device is used. The colors used to identify the products are mentioned below:

- Knotix PK is identified with a white and/or cyan blue (RGB 0-120-130) color,
- Knotix SPD PK is identified with a white and/or green (RGB 39-159-0) color,


- SpinFit PK features a white, orange (RGB 255-130-0) and/or green (39-159-0) handle, and
  - QuickFit PK comes in white, green (RGB 39-159-0) and/or orange (255-130-0).
  - Capillary: manufactured with stainless steel, as per the AISI 304 standard.
  - Rod: manufactured with stainless steel, as per the AISI 304 standard.
- Check the packaging labels for the anchor materials.

### 13. MORE INFORMATION

To learn more about this product, visit our website: [www.fergusmedical.com](http://www.fergusmedical.com).

### 14. SYMBOLS USED IN LABELS

	<b>CATALOG NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>SERIAL NUMBER</b>
	<b>PRECAUTION</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>DATE OF MANUFACTURE</b>
	<b>STERILE STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>MANUFACTURER</b>
	<b>SEE INSTRUCTIONS FOR USE</b>
	<b>DO NOT USE IF PACKAGING IS DAMAGED</b>
	<b>DO NOT RESTERILIZE</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>KEEP DRY</b>
	<b>MR-SAFE</b>

Fergus

# Fergus

**PROMEDON**

Av. Gral. Manuel Savio s/n  
Lote 3 • Manzana 3 • (X5123XAD)  
Pque Ind. Ferreyra, Córdoba, Argentina

[www.fergusmedical.com](http://www.fergusmedical.com)