Knotix Ti

Knotix Ti 2.8

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ESPAÑOL

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PORTUGUÊS

Instructions for use

ENGLISH

M-351-00204 (Rev. 02) / 27.04.2022

Fergus

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1. DESCRIPTION OF THE DEVICE

Knotix Ti and Knotix Ti 2.8 suture anchors are titanium devices designed with threads and eyelets for the passage of sutures. These products are supplied sterile and preloaded with non-absorbable sutures onto a disposable insertion device to facilitate implant placement.

Knotix Ti and Knotix Ti 2.8 titanium suture anchors are designed to provide fixation of soft tissue to bone.

The Knotix Ti model has an inner eyelet for the passage of CordFit #2 suture threads (USP standard). Two suture threads, either in blue, white and/or white/blue, can be threaded through the eyelet. The four ends of these two sutures are then secured over the suture retention area (see figure) located in the handle of the insertion device. This anchor is available in diameters of 5.5 mm or 6.5 mm and a length of 16.3 mm, which is included in the product label.

Figure 1: Knotix Ti anchor (illustrative picture).

1: Suture retention area.

2: Titanium anchor.

The Knotix Ti 2.8 model includes one CordFit #2 suture (USP standard), either in blue, white or white/blue, which is threaded through the anchor's outer eyelet. The ends of this suture are then secured over the retention area in the handle of the insertion device. This anchor is available in a diameter of 2.8 mm and a length of 11.6 mm, which is included in the product label.

Figure 2: Knotix Ti 2.8 anchor (illustrative picture).

1: Suture retention area.

2: Titanium anchor.

2. HOW SUPPLIED

Contents:

- **Knotix Ti Model:** implantable titanium anchor with two CordFit #2 sutures and an insertion device.
- **Knotix Ti 2.8 Model:** implantable titanium anchor with one CordFit #2 suture and an insertion device.

Find the model, anchor size, materials, and suture and insertion device color on each insertion device label.

All components are supplied sterile and ready-to-use.

3. INDICATIONS FOR USE

Knotix Ti and Knotix Ti 2.8 titanium suture anchors are intended to provide fixation of soft tissue to bone within the body joints. They are designed for the following surgical interventions:

- **Shoulder joint:**

- Rotator cuff repair
- SLAP lesion repair
- Bankart lesion repair

- **Ankle joint:**

- Lateral instability
- Medial instability

The adequate implant model and size must be determined at the surgeon's professional discretion based on the specific indication and the patient's medical history. However, for the shoulder joint, the Knotix Ti suture anchor (in either of its two diameters) is recommended for the repair of the rotator cuff (because of its size, it is more appropriate for the fixation of these tendons).

The Knotix Ti 2.8 suture anchor is recommended for the repair of SLAP and Bankart lesions (being smaller facilitates its implantation in the glenoid fossa).

For the repair of lateral and medial ankle instability, the 5.5-mm Knotix Ti and Knotix Ti 2.8 models can be used, depending on the physician's preference and the patient's condition.

The anchor provides an anchoring point for the suture on the bone surface, whereas sutures allow for tendon fixation to the anchoring point.

The insertion device allows for anchor implantation and suture placement in the right point.

4. CONTRAINDICATIONS

The device shall not be implanted in the following cases:

- Insufficient bone quality or quantity for tissue fixation. The physician should assess bone quality carefully in patients with lack of bone maturity before starting the orthopedic surgical procedure. The growth plate should not be deviated, interrupted or altered.
- Poor quality or pathological conditions of the soft tissue to be repaired, which could impair secure fixation by suture.
- Limited blood flow and/or previous infections which may delay healing.
- Conditions which tend to limit the patient's ability or will to follow doctor's instructions on postoperative care, which could include significant restrictions of their physical activity.

- Surgical interventions not listed in section “Indications for Use.”
- Inadequate neuromuscular status.

5. WARNINGS

- Before using this product, read the instructions for use fully and thoroughly.
- The device must be used only by a medical professional trained in the surgical techniques indicated for the application of these implants.
- The sale of this device is restricted, by law, to physicians or to anyone authorized by a physician.
- DO NOT use the product if the packaging is open or damaged, as product sterilization may have been compromised.
- Products are designed for a SINGLE use. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could damage their performance and increase the risk of inadequate sterilization and cross contamination.
- Maximum precautions must be taken to avoid contamination. Operating room conditions must meet hospital, administrative and/or local government requirements.
- Do not use the device after the expiration date stated on the label. The performance, safety and/or sterility of the device cannot be guaranteed after the expiration date.
- Ensure that the product is disposed of following local regulations and be aware of the environmental contamination risks to the staff and the patients.
- The product must be stored in its original packaging.
- The patient’s sensitivity to device materials must be considered before implantation.
- All surgical instruments used in the intervention meant for the application of this implant must be sterile and fit for use in aseptic environments.
- Product presentation should not be altered without following the surgical technique, as this could impair the implant’s performance.
- If the implant disengages from the insertion device, DO NOT REASSEMBLE, since this could damage its performance and increase the risk of inadequate sterilization.
- After the intervention and until healing is complete, the fixation provided by this device should be protected. The postoperative therapeutic guidelines prescribed by the physician must be followed thoroughly in order to prevent overstresses on the device and/or the bone.
- The patient will receive detailed instructions on the use and limitations of this device.
- If the explantation of the device is decided, the risks associated with the possibility of the patient undergoing a second surgical intervention must be

considered. The device removal must be combined with the appropriate surgical and postoperative approach.

- Open devices and implants should be discarded, even if they have not been used or handled, since this could increase the risk of inadequate sterilization and cross contamination.
- Biohazardous waste such as explanted devices, sutures, needles and contaminated surgical material must be disposed of safely and following the policy of the medical facility.
- The company is not responsible for failures or damage caused when these devices are combined with products not recommended by the company.
- The company shall not be liable for any effect or consequence that may result from not following specific surgical techniques or instructions, such as a mistaken indication, an inadequate use of a surgical technique, asepsis issues, etc., which could affect patient safety and device performance.

6. PRECAUTIONS

- The adequate implant model, quantity and size should be determined at the surgeon's professional discretion based on the specific indication and the patient's medical history.
- Surgeons are recommended to check the specific surgical technique for the product before proceeding with the operation.
- For Knotix Ti 2.8, when the bone presents a high resistance to insertion, drilling a hole on the bone surface might be necessary to avoid anchor damage.
- Do not remove the suture until the anchor has been completely inserted into the desired site, since this could impair implantation.
- Sutures can break if handled or managed with sharp or cutting instruments.
- Consult the company for the recommended placement/implantation systems to implant the device correctly.
- The decision of removing the implant shall be made considering the potential risk of a second surgical intervention for the patient. After implant removal, proper postoperative control is required.
- Prior to use, inspect the device in order to ensure that it is not damaged. Do not use it if it is damaged.
- Preoperative and postoperative procedures, including the knowledge of the surgical technique and the adequate selection and fixation of the implant, are important for the successful use of this implant.
- If sensitivity to materials is suspected, the appropriate tests must be carried out in order to rule out sensitivity before implantation.

7. SUPPLY AND STORAGE

- Knotix Ti and Knotix Ti 2.8 suture anchors are supplied sterile (using ethylene oxide) and pyrogen-free.
- In order to preserve its quality and efficacy, the product is supplied in a double Tyvek pouch. After sterilization, the product is placed in a cardboard box (outer packaging) to protect it from light and heat. The product is identified with labels. To avoid damage to the Tyvek pouches, a plastic protector is used, containing the insertion device and covering the implant.
- Keep the product away from any moisture or heat sources.
- The product is recommended to be kept at room temperature.
- DO NOT USE THE PRODUCT IF THE LABEL IS DAMAGED OR THE PACKAGING IS OPEN OR ALTERED.
- DO NOT USE AFTER THE EXPIRATION DATE STATED ON THE PACKAGING.
- The implant shelf life has been determined to be 3 years.

8. MAGNETIC RESONANCE (MR) ENVIRONMENT

1. MR conditional

Knotix Ti and Knotix Ti 2.8 titanium suture anchors by Fergus, made of Ti-6Al-4V ELI titanium alloy as per ASTM F136, feature a low magnetic susceptibility and electric conductivity. Therefore, there is a low risk of displacements due to induced magnetic forces or torsion, and of heating induced by radiofrequency or imaging devices.

Patients who have this device implanted can be subjected to magnetic resonance studies in a safe manner under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only.
- Magnetic field with a maximum spatial gradient of 3000 Gauss/cm.
- Maximum reported by the Magnetic Resonance System, specific absorption rate (SAR) averaged all over the body at 1 W/kg for 15 minutes of exploration in the normal mode of the MR system.

Image quality may be compromised in the area or the surroundings of the implant. An optimization of the magnetic resonance parameters is necessary in order to compensate for the presence of the device.

9. POTENTIAL COMPLICATIONS

Any surgical intervention involves a certain level of risk. Some complications that may occur are:

- Deep or superficial infections
- Atypical reactions to foreign bodies

- Pain at the incision or intervention site
- Implant loosening or detachment
- Bone injuries
- Incomplete or inadequate lesion repair or lesion recurrence

10. INFORMATION FOR PATIENTS

Surgeons are responsible for informing patients or their representatives, prior to the surgery, of the potential complications associated to the device implantation (see section "Potential complications").

The company and their distributors delegate to the surgeon the responsibility to inform the patient of the advantages and possible risks related to the implantation and use of the device.

Postoperative care is important. The patient will be instructed by the physician on the limitations of the implant and will be warned about the care needed for healing.

The patient must be warned that, if the explantation of the device is decided, the risks associated with the possibility of the patient undergoing a second surgical procedure must be considered. The device removal must be combined with the appropriate postoperative approach.

The performance of additional surgery does not guarantee the total repair of complications.

11. SURGICAL TECHNIQUE

The surgeon must choose the appropriate approach depending on the repair needed. For a foot or ankle repair, an open or arthroscopic approach may be used. For shoulder surgeries, the arthroscopic approach is recommended.

Additionally, the surgeon must choose the number of implants and the appropriate fixation layout, depending on the lesion type, the patient's medical history and other criteria deemed to be important for the successful use of the product in the surgical intervention.

The following instructions for surgical techniques only refer to the device implantation; they do not include other procedures that are essential to perform a surgical intervention, such as patient preparation, access portals, surgery completion, among others.

Knotix Ti

The Knotix Ti anchor is recommended to perform rotator cuff repairs in the shoulder and lateral and medial ankle instability repairs.

Steps for the right implantation of the Knotix Ti anchor:

1) Position the anchor in the desired implantation site.

2) Rotate the (insertion device) handle clockwise to enable anchor implantation in bone tissue by threading. While turning the handle, a slight forward pressure (in the proximal and distal direction) is recommended to make implantation easier.

ATTENTION: It is important to ensure that the anchor has been fully inserted to avoid possible damages to surrounding tissues and/or implant detachments. To this end, threading should continue until the laser line located in the distal end of the insertion device shaft has fully entered the hole where the anchor is placed.

3) Release sutures from the handle retention area and remove the insertion device.

4) Use a suture passer to pass sutures through the tissue, tie them and cut off the excess suture.

ATTENTION: If an open approach is used, the suture must be threaded through a needle to pass through the tissue, which allows for tissue fixation to the anchoring point. Tie surgical knots to end the fixation.

Note: The needle is not included in the kit.

After use, dispose of the product, sutures, needles, contaminated disposable surgical material and packaging, in accordance with hospital, administrative, or local government procedures.

Knotix Ti 2.8

The Knotix Ti 2.8 anchor is recommended for glenoid labrum repairs to correct shoulder instability and for lateral and medial ankle instability repairs.

Steps for the right implantation of the Knotix Ti 2.8 anchor:

1) Identify the anchor insertion site. Only in the case of shoulder instability repair, position an implantation guide of at least 3.2 mm of inner diameter in the insertion point.

2) Insert the anchor into the guide and rotate the (insertion device) handle clockwise to enable anchor implantation in bone tissue by threading. While turning the handle, a slight forward pressure (in the proximal and distal direction) is recommended to make implantation easier.

ATTENTION: It is important to ensure that the anchor has been fully inserted to avoid possible damages to surrounding tissues and/or implant detachments. To this end, threading should continue until the laser line located in the distal end of the insertion device shaft has fully entered the hole where the anchor is placed.

Note: If an open approach is used, no guide is required.

3) Release sutures from the handle retention area and remove the insertion device.

4) Use a suture passer to pass sutures through the tissue, tie them and cut off the excess suture.

ATTENTION: If an open approach is used, the suture must be threaded through a needle to pass through the tissue, which allows for tissue fixation to the anchoring point. Tie surgical knots to end the fixation.

Note: The needle is not included in the kit.

After use, dispose of the product, sutures, needles, contaminated disposable surgical material and packaging, in accordance with hospital, administrative, or local government procedures.

12. MATERIAL SPECIFICATIONS

Knotix Ti and Knotix Ti 2.8 suture anchors are made of Ti-6Al-4V ELI titanium alloy as per ASTM F136.

CordFit #2 sutures (USP standard) are made of ultra-high-molecular-weight polyethylene (UHMWPE). Sutures can be blue, white or white/blue in color. The color additive 73.1015 chromium-cobalt-aluminum oxide is the pigment used to dye blue sutures and the blue portion of the white/blue sutures; no dye is used for white sutures.

The insertion device consists of two parts:

- **Handle:** It is made of medical grade ABS and it has different colors to identify the anchor model for which this insertion device is used. A white and/or cyan blue (RGB #00828c) device is used for a Knotix Ti suture anchor, whereas a white and/or green (RGB #a7d500) insertion device is used for a Knotix Ti 2.8 suture anchor.
- **Shaft:** It is made of stainless steel as per AISI 304.

See the packaging label to learn more about the anchor materials.

13. MORE INFORMATION

If you need more information about this product, visit our web page, www.fergusmedical.com.

14. SYMBOLS USED IN THE LABELS

	CATALOGUE NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	PRECAUTION
	EXPIRATION DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	SEE INSTRUCTIONS FOR USE
	DO NOT USE IF THE PACKAGING IS DAMAGED
	DO NOT RESTERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	MR CONDITIONAL

Fergus

Fergus

PROMEDON SA

Av. Gral. Manuel Savio s/n

Lote 3 • Manzana 3 • (X5123XAD)

Pque. Ind. Ferreyra, Córdoba, Argentina

www.fergusmedical.com