**Full** Blade  
**Master** Blade  
**Total** Blade  
**Cyl** Burr  
**Ball** Burr

---

Puntas de Shaver

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Instructions for use  
ENGLISH

M-350-00845 (Rev. 02) / 1.11.2022

Fergus

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## ENGLISH

### 1. DEVICE DESCRIPTION

Fergus shaver blades and burrs are specifically designed with different shaver/burr configurations and cutting windows for soft tissue, cartilage and bone resection during arthroscopic surgical procedures.

The Fergus shaver blades and burrs include the Full Blade, Master Blade, Total Blade, Cyl Burr and Ball Burr models.

The Full Blade model is designed for the meniscal, cartilage, and synovial membrane and plicae debridement.

The Master Blade model is designed for soft tissue debridement in the knee and the shoulder.

The Total Blade model is designed for articular cartilage debridement as well as for the resection of the periosteum of the acromion.

The Cyl Burr model is designed for swift and aggressive resection of the femoral bone in intercondylar notch remodeling, or of the distal clavicle, and for subacromial decompression procedures.

The Ball Burr model is designed to be used for bone resection of the intercondylar notch and subacromial decompression.

## **2. PRODUCT PRESENTATION**

The product comes in a box containing 5 units of the same model, each of them packed in a blister.

Check the model and the diameter of the shaver blade/burr in each of the individual device labels.

## **3. INDICATIONS FOR USE**

The Fergus shaver blades and burrs are designed for soft and bone tissue resection during arthroscopic surgical procedures.

## **4. CONTRAINDICATIONS**

- Insufficient bone quantity or quality: the physician shall carefully assess the bone quality of the patients with bone immaturity before starting the orthopedic surgical procedure. The growth cartilage must not be deviated, interrupted or altered.
- Limited blood irrigation and/or previous infection which may delay consolidation.

## **5. WARNINGS**

- Read the instructions for use carefully and completely before using the product.
- The device must be used only by a trained medical professional.
- The use of this device is restricted to physicians or people authorized by a physician.
- DO NOT use the device if the packaging is open or damaged since product sterility will be lost.
- Shaver blades and burrs have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their efficiency and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative or governmental requirements.
- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- The product must be stored in the original packaging.
- Bio-hazardous waste and contaminated surgical

materials must be disposed of safely and as per the policies of the site where the surgical procedure is performed.

- The manufacturer shall not be liable for any effects or consequences which may arise from deviating from these specific instructions.

## **6. PRECAUTIONS**

- Surgeons must follow their professional judgement to determine the appropriate device model to use for the surgery.
- Device components must not be interchanged.
- Do not allow any rotating parts of the device to come into contact with metal parts, since this could cause irreversible damage to both objects. Contact between metallic parts could change the sharpness of the cutting edge of the blade and/or cause the tip to break inside the body.
- Do not apply excessive force to the device.
- Do not use the device without irrigation flow. Otherwise, you will increase the temperature of the device and, therefore, cause burns to the patient. Do not handle the cutting terminal end.

## **7. SUPPLY AND STORAGE**

- The devices are supplied sterile. They are sterilized using ethylene oxide. In order to maintain product quality and effectiveness they are marketed in a box with the appropriate labels. Each unit is also contained in a blister.
- Keep the products away from any sources of humidity.
- Room temperature storage is advised for this product.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- The device has a 3-year shelf life.

## **8. POSSIBLE COMPLICATIONS**

- Additional damage due to applying excessive force while using the device
- Infection, both deep and superficial

## **9. ASSEMBLY**

Open the packaging, remove the sterile product carefully, and place it in a sterile area to check its integrity.

Mount the shaver blade/burr to the handpiece following these instructions:

### **a. Coupling:**

To place the shaver blade/burr in the handpiece, it must be inserted so that the locking structure is aligned with the handpiece cavity. If it does not become locked in place successfully, remove the shaver blade/burr, turn its internal shaft to align

it with the engine output shaft and place it again, making sure it locks securely.

**b. Decoupling:**

To remove the device from the handpiece, press the button in the handpiece to decouple the blade/burr.

**10. MATERIAL SPECIFICATIONS**

Fergus shaver blades and burrs are made with AISI 304 and AISI 455 stainless steel, polycarbonate, 6351 aluminum and polytetrafluoroethylene.

**11. LEARN MORE**

To learn more about this product, visit our website: [www.fergusmedical.com](http://www.fergusmedical.com)

**12. SYMBOLS USED IN LABELS**

	<b>CATALOG NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>SERIAL NUMBER</b>
	<b>PRECAUTIONS</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>MANUFACTURING DATE</b>
	<b>STERILE STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>MANUFACTURER</b>
	<b>CHECK INSTRUCTIONS FOR USE</b>
	<b>DO NOT USE IF PACKAGING IS DAMAGED</b>
	<b>DO NOT RE-STERILIZE</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>KEEP DRY</b>

# Fergus

**PROMEDON SA**

Av. Gral. Salvo s/n

Lote 3 • Manzana 3 • (X5123XAD)

Pque Ind. Ferreyra, Córdoba, Argentina

[www.fergusmedical.com](http://www.fergusmedical.com)