# Spirom Ti

# Spirom Ti BTB

# Spirom Ti FT

---

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## ENGLISH

### 1. DESCRIPTION OF THE DEVICE

The Fergus titanium interference screws include the **Spirom Ti** model, the **Spirom Ti BTB** model, and the **Spirom Ti FT** model.

The **Spirom Ti** model is designed to fix semitendinosus muscle tendon grafts, gracilis muscle tendon grafts, sartorius muscle tendon grafts, and quadriceps muscle tendon grafts.

The **Spirom Ti BTB** model is designed to fix bone-patellar tendon-bone grafts.

The **Spirom Ti FT**, with a full thread, is designed to fix semitendinosus muscle tendon grafts, gracilis muscle tendon grafts, sartorius muscle tendon grafts, and quadriceps muscle tendon grafts.

Interference screws are threaded and cannulated.

### 2. HOW SUPPLIED

Contents:

1 unit Titanium interference screw.

Find the interference screw model, diameter and length on the individual device labels.

### 3. INDICATIONS FOR USE

These devices are designed to provide interference fixation in anterior and posterior knee cruciate ligament reconstruction.

### 4. CONTRAINDICATIONS

- Insufficient bone quality or quantity for attachment. The physician must evaluate carefully the bone quality of patients with bone immaturity before starting with the orthopedic surgical procedure, and must avoid deviating, interrupting or altering the growth cartilage.
- Limited blood irrigation and/or previous infection which may delay consolidation.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions on their physical activity.

## 5. WARNINGS

- Before using this product, read the instructions for use fully and thoroughly.
- The device must be used only by a medical professional trained in the surgical techniques indicated for the application of these implants.
- The law restricts the use of this device to physicians or people authorized by a physician.
- DO NOT use the device if the packaging is open or damaged since product sterility will be lost.
- The device must not be used for surgical procedures that are not included in the "Indications for Use" section.
- Internal fixation devices are designed for a SINGLE use. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could damage their performance and increase the risk of inadequate sterilization and cross contamination.
- Maximum precautions must be taken to avoid contamination. Operating room conditions must meet hospital, administrative and/or local government requirements.
- Do not use the device after the expiration date stated on the label. The performance, safety and/or sterility of the device cannot be guaranteed after the expiration date.
- The product must be stored in its original packaging.
- The patient's sensitivity to device materials must be considered before implantation. See section "Potential complications."
- The screw may break or fixation may be lost if the insertion site is not properly prepared.
- All the surgical instruments used in the surgical intervention meant for the application of this implant must be sterile and fit for use in aseptic environments.
- Before performing the insertion, inspect and confirm that the interference screw fits completely on the screwdriver. If the interference screw does not fit correctly on the screwdriver, the screw may break or get damaged.
- DO NOT use a different screwdriver from the one stated in the "Surgical Technique" section.
- DO NOT use a different Nitinol guide from the one stated in the "Surgical Technique" section.
- After the operation and until final consolidation, the fixation provided by the device must be considered as temporary, and therefore it might not withstand weight loads or other stresses without assistance. The fixation provided by the device must be protected. The postoperative therapeutic guidelines prescribed by the physician must be followed thoroughly in order to prevent overstresses on the device and/or the bone.
- If the explantation of the device is decided, the risks associated with the possibility of the patient undergoing a second surgical procedure must be considered. The device removal must be combined with the appropriate postoperative approach.
- In the case of metal implants, the use of screw extraction instruments may be required if the devices have

been implanted for a prolonged period of time.

- Ensure that the product is disposed of following local regulations and be aware of the environmental contamination risks to the staff and the patients.
- Biohazardous waste such as explanted devices, needles, sharp objects and contaminated surgical material must be disposed of in a safe manner and following the policy of the hospital where the surgical intervention is performed.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not authorized/recommended by the company.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from the specific surgical techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, or problems with asepsis.

## 6. PRECAUTIONS

- The adequate screw model and size must be determined at the surgeon's professional discretion based on the specific indication, the preferred surgical technique, and the patient's medical history.
- Surgeons are recommended to check the specific surgical technique for the device before proceeding with the operation.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.
- The preoperative and postoperative procedures, including an adequate selection of surgical techniques and the correct selection and fixation of the implant, are important considerations when using this device.
- Before using, visually inspect the device for damage. Do not use if damaged.
- Inspect the surgical screwdriver before inserting the screw, in order to ensure there are no scratches, scrapes or marks. Make sure that the end is not bent in a way that could prevent the screw from fitting. Such defects increase the risks of breaking the screw and/or the screwdriver.
- The surgeon must ensure that there are no portions of the implant protruding after its insertion. Otherwise, soft tissue irritation or postoperative pain may occur.
- The knee position must be kept fixed. Its flexion angle must not be changed during screw insertion, since this may cause screw divergence or screwdriver malfunctioning.
- Excessive forces should be avoided on the screw, the bone and/or the insertion instruments.

## 7. SUPPLY AND STORAGE

- The devices are supplied sterile and pyrogen-free, with ethylene oxide as the sterilization means. In order to preserve the quality and efficacy of the products,

they are commercialized in a double Tyvek pouch. The implants must be used only if the factory packaging and labeling are intact.

- Keep the products away from any moisture sources.
- The products are recommended to be kept at room temperature.
- DO NOT USE THE PRODUCT IF THE PACKAGING IS OPEN OR DAMAGED.
- DO NOT USE AFTER THE EXPIRATION DATE STATED ON THE PACKAGING.
- The implant shelf life has been determined to be 3 years.

## **8. MAGNETIC RESONANCE (MR) ENVIRONMENTS**

- MR conditional:

The Fergus-brand titanium interference screw, manufactured with the Ti-6Al-4V ELI titanium alloy as per ASTM F136, exhibits low magnetic susceptibility and low electrical conductivity; therefore, the risks of displacement caused by magnetically-induced forces or torque, radiofrequency-induced heating and imaging artifacts decrease.

Patients who have this device implanted can be subjected to magnetic resonance studies in a safe manner under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only.
- Magnetic field with a maximum spatial gradient of 3000 Gauss/cm.
- Maximum reported by the Magnetic Resonance System, specific absorption rate (SAR) averaged all over the body at 1 W/kg for 15 minutes of exploration in the normal mode of the MR system.

Image quality may be compromised in the area or the surroundings of the implant. An optimization of the magnetic resonance parameters is necessary in order to compensate for the presence of the device.

## **9. POTENTIAL COMPLICATIONS**

- Infection, either deep or superficial.
- Mild inflammatory reaction.
- Allergic reaction.
- Reaction to foreign bodies.
- Arthrofibrosis.

## **10. PATIENT INFORMATION**

Surgeons are responsible for informing patients or their representatives in advance on the possible complications associated to the device implantation (see "Potential complications").

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated with the

implantation and use of the device.

Postoperative care is important. The patient will be given instructions on the limitations required by the implant and will be warned about weight bearing and overloading the device before safe bone consolidation. The patient must be warned that, if the explantation of the device is decided, the risks associated with the possibility of the patient undergoing a second surgical procedure must be considered. The device removal must be combined with the appropriate postoperative approach.

The performance of additional surgery does not guarantee the total repair of complications.

## 11. SURGICAL TECHNIQUE

The screw is designed to be used for graft fixation in knee cruciate ligaments reconstruction.

The following surgical technique instructions refer only to device implantation, not considering other procedures that are critical to perform the surgical procedure, such as patient preparation, working portals, procedure finalization, etc.

This implant is suitable to be fitted either on the tibial or femoral bone tunnel.

The femoral tunnel can be performed by means of the following surgical techniques:

- Transportal technique (through the anteromedial arthroscopic portal).
- Transtibial technique (through the graft tibial bone tunnel).
- Outside-in technique.

### **General Steps:**

1. Remove and prepare the graft using appropriate techniques and instruments.
2. Drill the bone tunnels with the appropriate instruments, taking into account the chosen surgical technique, the anatomic references and the graft dimensions.
3. Position the graft inside the joint in such a way that it fits the bone tunnels created.
4. The surgeon must choose the appropriate screw model and size taking into account the size of the bone tunnels created, the type and size of the graft, the indication for use of the product, the chosen surgical technique, and the patient's medical history.
5. Fit the interference screw on an appropriate surgical screwdriver with a 3.5 mm hexagonal tip and a minimum hexagon length of 5.5 mm, and position it on the bone tunnel end. A sterile guide wire of 2.0 mm must be used to mount the screw and the screwdriver on it, in order to guarantee a correct insertion into the tunnels.
6. Next, turn the screw clockwise in order to screw it, applying tension on the graft. The surgeon must screw until they consider that the fixation depth and tension between the graft and the bone walls are adequate.
7. Separate the screwdriver from the screw and remove the guide wire.

8. Repeat steps 5 and 6 for the remaining screw and bone tunnel.

### **Warnings**

- The graft diameter must be measured accurately to ensure that it corresponds to the bone tunnel diameter.
- If the implant is chosen to be fitted to the femur, care must be taken not to damage the graft while the screw is slid to its insertion position.
- If the screw fixation depth is not adequate, no appropriate interference fixation will be achieved.
- Do not use the screwdriver and screw assembly as a lever. The screwdriver and screw assembly must be aligned axially with the drill bit channel.

## **12. MATERIAL SPECIFICATIONS**

The interference screw is made of Ti-6Al-4V ELI titanium alloy as per ASTM F136. See the packaging label to learn more about the materials.

## **13. MORE INFORMATION**

If you need more information about this product, visit our web page [www.fergusmedical.com](http://www.fergusmedical.com)

## **14. SYMBOLS USED IN LABELS**

	<b>CATALOGUE NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>PRECAUTION</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>DATE OF MANUFACTURE</b>
	<b>STERILE. STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>MANUFACTURER</b>
	<b>DO NOT USE IF THE PACKAGING IS DAMAGED</b>
	<b>SEE INSTRUCTIONS FOR USE</b>
	<b>DO NOT RESTERILIZE</b>
	<b>SERIAL NUMBER</b>

**KEEP AWAY FROM SUNLIGHT**

**KEEP DRY**

**MR CONDITIONAL**

# Fergus

**PROMEDON SA**

Av. Gral. Manuel Savio s/n

Lote 3 • Manzana 3 • (X5123XAD)

Pque. Ind. Ferreyra, Córdoba, Argentina

[www.fergusmedical.com](http://www.fergusmedical.com)