# Spirom PK

---

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## 1. DESCRIPTION OF THE DEVICE

The Fergus-brand, **Spirom PK**-model PEEK interference screws are designed for the fixation of soft tissue grafts or bone-tendon-bone grafts to the tibia or the femur during anterior and posterior cruciate ligament reconstruction procedures in the knee. Interference screws are threaded and cannulated.

## 2. HOW SUPPLIED

Contents:

1 PEEK interference screw

Check the diameter and length of the PEEK interference screw in each of the individual device labels.

## 3. INDICATIONS FOR USE

These devices are designed to provide interference fixation for the reconstruction of the anterior or posterior cruciate ligament in the knee.

## 4. CONTRAINDICATIONS

- Insufficient bone quality or quantity for attachment. The physician must evaluate carefully the bone quality of patients with bone immaturity before starting with the orthopedic surgical procedure, and must avoid deviating, interrupting or altering the growth cartilage.
- Limited blood irrigation and/or previous infection which may delay consolidation.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions on their physical activity.

## 5. WARNINGS

- Read the instructions for use carefully and completely before use.

- The device must be used only by a medical professional trained in the surgical techniques indicated for placing these implants.
- The law restricts the use of this device to physicians or people authorized by a physician.
- DO NOT use the device if the packaging is open or damaged since product sterility will be lost.
- The device must not be used for surgical procedures that are not included in the "Indications for Use" section.
- Internal fixation devices have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their efficiency and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative or governmental requirements.
- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- The product must be stored in the original packaging.
- Patient sensitivity to device materials must be considered before implantation.
- The screw may break or fixation may be lost if the insertion site is not properly prepared.
- All surgical instruments used in the surgical procedure to install this implant must be sterile and suitable for use in aseptic environments.
- Check and confirm that the interference screw is completely held by the screwdriver before performing the insertion. If the interference screw is not properly held by the screwdriver, it may break or become damaged.
- DO NOT use a different screwdriver from the one stated in the "Surgical Technique" section.
- DO NOT use a different nitinol guide from the one stated in the "Surgical Technique" section.
- After the procedure and until complete consolidation is achieved, the fixation provided by the device must be considered as temporary. Therefore, it may not be able to resist weight loads or other efforts without assistance. The fixation provided by the device must be protected. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing the device and the bone.

- If the device is to be explanted, the risks of the possibility of a second surgical procedure on the patient must be considered. Device extraction must be accompanied by an appropriate post-operative approach.
- Make sure that the product is disposed as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- Bio-hazardous waste, such as explanted devices, needles, sharp objects and contaminated surgical materials, must be disposed of safely and as per the policies of the site where the surgical procedure is performed.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the company.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, which may affect patient safety and implant performance.

## 6. PRECAUTIONS

- Surgeons must follow their professional judgment to determine the appropriate screw size based on the specific indication, the preferred surgical technique and the history of the patient.
- Surgeons are advised to review the device-specific surgical technique before the procedure.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.
- Pre- and post-operative procedures, including the appropriate selection of surgical techniques and the selection and correct placement of the implant, are important considerations when using this device.
- Before using, visually inspect the device for damage. Do not use if damaged.
- Inspect the surgical screwdriver before inserting the screw in order to make sure that there are no scratches, scrapes or marks. Make sure that the tip is not bent in a way that could prevent it from fitting in the screw. These defects increase the risk of the screw and/or the screwdriver breaking.

- The surgeon must ensure that no parts of the implant stick out after implantation. If not, it could cause irritation in soft parts or postsurgical pain.
- Knee position must be kept immobilized. Knee flexion angle must not change during screw insertion, since it could cause screw divergence or screwdriver malfunctions.
- Avoid applying excessive force on the screw, the bone or the insertion instruments.

## **7. SUPPLY AND STORAGE**

- Devices are supplied sterile and pyrogen-free. They are sterilized using ethylene oxide. In order to maintain product quality and effectiveness, they are marketed in a double Tyvek pouch. Implants must be used only if the original packaging and labeling are not damaged.
- Keep the products away from any sources of humidity.
- Room temperature storage is advised for this product.
- **IF THE PACKAGING IS OPEN OR ALTERED, DO NOT USE THE PRODUCT.**
- **DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.**
- The implant has a 3-year shelf life.

## **8. MAGNETIC RESONANCE (MR) ENVIRONMENTS**

### **MR Safety**

The interference screws are manufactured solely with polyether ether ketone (PEEK); therefore, they are safe for MR.

## **9. POTENTIAL COMPLICATIONS**

- Infection, both deep and superficial.
- Mild inflammatory reactions.
- Allergic reactions.
- Reactions to foreign bodies.
- After surgery, ligament laxness can be affected, and therefore, the knee range of motion may decrease or increase.

## 10. PATIENT INFORMATION

The surgeon is responsible for informing the patient and his/her representatives before surgery about the possible complications related to the implantation of the device (see section "Potential complications").

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated with the implantation and use of the device.

Post-operative care is important. The patient will be instructed about the limitations required due to the implant and will be warned about weight unloading and device overloading before a safe bone consolidation.

The patient must be warned that, if explantation of the device is deemed necessary, the risks of the possibility of a second surgical procedure on the patient must be considered. Device extraction must be accompanied by an appropriate post-operative approach.

Performing additional surgeries does not guarantee a complete repair of the complications.

## 11. SURGICAL TECHNIQUE

The screw is designed to be used as fixation in knee cruciate ligament reconstruction procedures.

The following surgical technique instructions refer only to device implantation, not considering other procedures that are critical to perform the surgical procedure, such as patient preparation, working portals, procedure finalization, etc.

This implant is suitable to be placed either in the femoral or tibial bone tunnel.

The femoral tunnel may be performed by using the following surgical techniques:

- Transportal technique (through the anteromedial arthroscopic portal)
- Transtibial technique (through a tibial bone tunnel for the graft)
- Inside-outside technique

1. Remove and condition the graft by using the appropriate techniques and instruments.

2. Create the bone tunnels considering the chosen surgical technique, the anatomical preferences and the dimensions of the graft, and using appropriate instruments.

The diameter of the bone tunnel is recommended to be equivalent to that of the graft.

3. Position the graft inside the joint so that it matches the bone tunnels created.
  4. The surgeon must select the screw size considering these factors: dimensions of bone tunnels created, graft type and size, product indications for use, surgical technique selected and patient history.
  5. Place the interference screw on a hexalobular (6-point) sterile surgical screwdriver for PEEK interference screws, and place it at the end of the bone tunnel. The screw and the screwdriver must be mounted over a 1.1 mm sterile nitinol guide wire in order to allow for a proper placement into the tunnels.
- Note:** If a bone-tendon-bone graft is used, the surgeon may opt to perform the preparation (thread guide) of the tunnel/bone block before inserting the screw. The necessary instrument is a 7, 8, 9 or 10 mm male tap, depending on the dimensions of the screw to be inserted.
6. Insert the screw by turning the screwdriver clockwise. The surgeon must thread until reaching a depth and a fixation pressure between the graft and the bone wall that the surgeon deems appropriate.
  7. Detach the screwdriver from the screw and remove the guide wire.
  8. Repeat steps 5 and 6 for the remaining bone tunnel, if this fixation method is chosen.

**Warnings:**

- Graft diameter must be measured accurately so that it matches the diameter of the bone tunnel.
- If the implant is to be placed in the femur, carefully avoid damaging the graft while sliding the screw towards the insertion site.
- If the screw placement depth is not appropriate, a correct interferential fixation will not be achieved.
- Do not use the screwdriver-screw assembly as a lever. The screwdriver-screw assembly must be aligned axially with the drill bit flute.

## 12. MATERIAL SPECIFICATIONS

The interference device is manufactured with polyether ether ketone (PEEK), as per the ASTM F2026 standard.

Check the packaging labels for information about the materials.

### 13. MORE INFORMATION

To learn more about this product, visit our website:  
[www.fergusmedical.com](http://www.fergusmedical.com)

### 14. SYMBOLS USED IN LABELS

	<b>CATALOG NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>MANUFACTURE DATE</b>
	<b>STERILE STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>MANUFACTURER</b>
	<b>DO NOT USE IF PACKAGING IS DAMAGED</b>
	<b>CHECK THE INSTRUCTIONS FOR USE</b>
	<b>CAUTION</b>
	<b>DO NOT RESTERILIZE</b>
	<b>SERIAL NUMBER</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>KEEP DRY</b>
	<b>MR-SAFE</b>
	<b>DOUBLE STERILE BARRIER SYSTEM</b>

# Fergus

**PROMEDON S.A**

Av. Gral. Manuel Savio s/n

Lote 3 • Manzana 3 • (X5123XAD)

Pque Ind. Ferreyra, Córdoba, Argentina

[www.fergusmedical.com](http://www.fergusmedical.com)