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1. DEVICE DESCRIPTION

Fergus Monopolar Electrodes are made of a terminal with a live tip followed by a handle with two buttons to select the operation mode, and finally, a cable with a 3-pin connecting terminal. There are different models that share these characteristics but use different tips.

The **FlatGrid** model features a grid tip and comes in two versions: rigid and flexible.

The **HookedShape** model has a hook-shaped tip and comes in two versions: rigid and flexible.

The **ProCurved** model includes a convex tip.

The **InCurved** model comes with a concave tip.

2. PRODUCT PRESENTATION

The product comes in a box containing (1) unit of the requested model. Each unit is packed in a blister.

Check the Monopolar Electrode model in the individual device label.

3. INDICATIONS FOR USE

Fergus Monopolar Electrodes are designed for resection, ablation, tissue removal, blood vessel hemostasis and soft tissue coagulation. They are used in knee, shoulder, ankle, wrist and elbow procedures.

4. CONTRAINDICATIONS

- Do not use the Monopolar Electrodes in any surgical procedures other than arthroscopic procedures.
- Do not use in patients where arthroscopic procedures are contraindicated for any reason.
- They must also not be used during surgical procedures without connecting a return electrode to the patient.
- Do not use in patients with pacemakers or other electronic implants.
- Do not use in arthroscopic procedures where non-conductive liquid media are used.

5. WARNINGS

- Read the instructions for use carefully and completely before use.
 - The product must be operated at temperatures between 10°C and 40°C, with a relative humidity between 30% and 75%, without condensation, and with an atmospheric pressure between 700 hPa to 1060 hPa.
 - The use of this device is restricted to physicians or people authorized by a physician.
 - Fergus Monopolar Electrodes have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their efficiency and increase the risk of inappropriate sterilization and cross-contamination.
 - The tip of the electrode must be constantly surrounded with conductive liquid. Do not use non-conductive irrigation liquids (sterilized water, air, gas, glycine, purisol, sorbitol-manitol, etc.).
 - DO NOT use the device if the packaging is open or damaged since product sterility will be lost.
 - Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
 - The product must be stored in the original packaging.
 - Do not bend the device, since this could cause the rupture and deterioration of the isolating material.
 - Do not use the electrode with a metallic cannula. The use of non-conductive cannulas is advised.
 - Avoid contact between the device and the arthroscope, since it could cause damage to the patient and/or devices.
 - Do not use when flammable materials are present. Do not place devices that carry an electrical current near or in contact with flammable materials (gauze, surgical drapes, etc.). This could cause a fire.
 - Do not touch the tip of the electrode if it is energized.
 - Do not insert or remove the live tip from the place of the surgery while the device is activated.
 - Excessive monopolar electrode wear may result from vigorous use against bone tissue.
- The tip of the device may remain hot and cause burns even after deactivating the device.
- Bio-hazardous waste and contaminated surgical materials must be disposed of safely and as per the policies of the site where the surgical procedure is performed.
 - The device should not be used in surgical procedures beyond those stated in the "Indications for Use" section.
 - The manufacturer shall not be held responsible for faults or damages caused when combining these devices with products not recommended by the company.
 - The manufacturer shall not be liable for any effects or consequences which may arise from deviating from these specific instructions.

6. PRECAUTIONS

- Inspect the device packaging. If it is damaged, do not use the device.
- Inspect the device visually before use. If there are visual defects or isolating material deterioration, do not use this device and replace it with another.
- Check the return plate for proper operation and ensure that the placement and the surface in contact with the patient are correct.
- Set the operation mode carefully, using the least amount of power possible and the least contact time needed for a satisfactory procedure.
- Using high power –over 90 W– and long exposure time values may cause damage to the isolating material, fuse the tip and even cause injuries and burns to the patient.

7. SUPPLY AND STORAGE

- The devices are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide. In order to maintain product quality and effectiveness it is supplied in a blister. After the sterilization, the product is protected during transportation and storage by placing it inside a cardboard box. Labels may be added for identification.
- Keep the products away from any sources of heat and humidity.
- We recommend storing the product at temperatures between 10°C and 40°C, with a relative humidity between 30% and 75%, without condensation, and with an atmospheric pressure between 700 hPa to 1060 hPa.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- Product shelf life is 3 years.

8. POSSIBLE COMPLICATIONS

- Burns
- Lesions to adjacent structures and tissue
- Interference with pacemakers and implantable electronic devices
- Conversion to open surgery
- Severe postoperative pain
- Persistent bleeding
- Infection
- Surgery delay

9. MODE OF USE

- Prepare the patient using electrosurgical procedure techniques.
- Carefully open the outside packaging and place the sterile device on the sterile field.
- Before using, visually inspect the device for damage. Do not use it if damaged.
- Insert the monopolar connector properly into the electrosurgical generator.

WARNING: During surgery, never leave the monopolar electrode outside the field of vision to prevent accidental tissue damage.

- Turn on the electrosurgical generator.
- Check for proper setting of electrosurgical controls before surgery: power setting (cutting, coagulation) compared with the maximum values declared for the monopolar electrode.

The power settings are listed below for the specific modes of use and technical characteristics (These values may be changed based on professional judgment, but the limit should be 90 W for both modes):

- Power Settings

SETTING/MODE	CUTTING	COAGULATION
RECOMMENDED POWER	30 W maximum	35 W maximum
MAXIMUM NOMINAL VOLTAGE	348 V	744 V

- Technical Specifications
 - Maximum peak voltage: 90 kV
 - Maximum output current: 1 A
 - Frequency limits: 0-1.0 MHz
 - Resistance: 300 Ω
 - Output power: 50-80 W
- Perform the procedure with the Monopolar Electrode controlling the speed, duration, closeness to tissue, and the setting of the buttons for cutting and coagulation as per the objectives and the specific situation in the surgical context.
- To use the cutting function: press the yellow button on top of the device handle.
- To use the coagulation function: press the blue button on top of the device handle.
- To remove residual necrotized tissue (scabs) from the tip of the Monopolar Electrode: remove the device from the patient, rub the tip against an abrasive cloth and then against a piece of wet gauze.
- Turn off the electrosurgical unit and remove the connector from the electrosurgical generator.
- Discard the monopolar electrode by following the protocols set by the hospital.

- Regardless of the type of surgery, always keep the tip of the device submerged in a conductive irrigation solution (saline solution, Ringer's lactate, etc.). An appropriate irrigation flow must be maintained to prevent the liquid from possibly overheating; if not, the patient may suffer heat-related injuries.




10. MATERIAL SPECIFICATIONS

Fergus Monopolar Electrodes are made up of an stainless steel shaft and tip. The tip features a protection made of ceramic material, and the remainder of the tip and the shaft are covered with polyamide paint for its protection. The handle of the Monopolar Electrode is made of acrylonitrile butadiene styrene.

11. LEARN MORE

To learn more about this product, visit our website:
www.fergusmedical.com

12. SYMBOLS USED IN LABELS

	CATALOG NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	PRECAUTION
	EXPIRATION DATE
	DO NOT REUSE
	MANUFACTURING DATE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE SINGLE STERILE BARRIER SYSTEM
	MANUFACTURER
	CHECK INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	DO NOT RE-STERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	STORAGE AND TRANSPORTATION TEMPERATURE LIMITS: 10°C TO 40°C
	STORAGE AND TRANSPORTATION RELATIVE HUMIDITY LIMITS: 30% TO 75%
	STORAGE AND TRANSPORTATION ATMOSPHERIC PRESSURE LIMITS: 700 hPa TO 1060 hPa

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PROMEDON S.A.