Cánulas Rígidas

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Instructions for use
ENGLISH

Fergus

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1. DEVICE DESCRIPTION

Fergus cannulas are supplied together with a disposable obturator that allows for their insertion in the joint space. The cannulas are threaded, to prevent movement when the instrument is removed, and transparent, which enables the user to see the instruments and the implantable material introduced through them.

They also feature a lateral portal with a valve for irrigation to go in or liquid to come out.

Fergus cannulas are available in different diameters and lengths, with color coding to set them apart.

2. PRODUCT PRESENTATION

The product comes in boxes containing five (5) or six (6) units. Each of the units is contained in an individual blister. Check the Cannula model in the individual device label.

3. INDICATIONS FOR USE

The Fergus rigid clear cannulas are indicated to be used as entry and exit portals in arthroscopic surgeries. During arthroscopic procedures, they allow for direct visualization and passage of instruments and sutures, and prevent water spillages. Additionally, their lateral channel can be used for irrigation entry, passive exit or controlled exit.

4. CONTRAINDICATIONS

- Sensitivity and/or adverse reactions to the materials that make up the device. If an adverse reaction is suspected, perform sensitization tests.

5. WARNINGS

- Read the instructions for use carefully and completely before use.

- The use of this device is restricted to physicians or people authorized by a physician.
- DO NOT use the device if the packaging is open or damaged since product sterility will be lost.
- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- The product must be stored in the original packaging.
- This is a single-use device. Therefore, its cleaning and resterilization have not been assessed, and, for this reason, its reuse may negatively affect product safety and effectiveness.
- Bio-hazardous waste and contaminated surgical materials must be disposed of safely and as per the policies of the site where the surgical procedure is performed.
- The device should not be used in surgical procedures beyond those stated in the “Indications for Use” section.
- The use of reusable obturators from a non-Fergus brand is not validated. Therefore, this could compromise the safety and effectiveness of the device.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the company.
- The manufacturer shall not be liable for any effects or consequences which may arise from deviating from these specific instructions.

6. PRECAUTIONS

- Inspect the device packaging. If it is damaged, do not use the device.
- In order to preserve seal integrity and prevent leaks, do not use devices with a diameter greater than that of the cannula.
- Applying flexion forces on the device while inserted in the joint may cause it to break or warp.
- Ensure a full insertion of the device into the joint by making sure it is completely visible inside the joint cavity and that it has penetrated through the appropriate layers of tissue. This reduces the probability of liquid leaks due to the tip being inserted in soft tissue.

7. SUPPLY AND STORAGE

- The devices are supplied sterile. They are sterilized using ethylene oxide. In order to maintain product quality and effectiveness, they are marketed in a box with the appropriate labels. Each unit is also contained in a blister.
- Keep the products away from any source of humidity.
- Room temperature storage (between 10 °C and 35 °C) is advised for this product.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- The device has a 3-year shelf life.

8. POSSIBLE COMPLICATIONS

- Infection
- Adverse reaction to the materials that make up the device

9. MODE OF USE

Insertion

- Perform an incision with a scalpel at the point of the skin in the anatomical region where you want to enter.
- Remove the elements from the sterile packaging using aseptic techniques.
- Insert the disposable obturator into the cannula checking that the obturator protrusions are aligned with the cannula grooves. Note: if a reusable obturator is used, consider the diameter and the length of the cannula for its selection. (Fergus reusable obturators are the same color as their corresponding cannula.) Insert a metal guide to define the placement site in the joint more precisely. Afterwards, mount over this guide the appropriate obturator previously loaded with the cannula.
- Insert the cannula-obturator assembly into the incision by applying continuous pressure downwards and threading clockwise.
- Through arthroscopic visualization, check that the cannula has been completely inserted.
- Remove the obturator while maintaining the position of the cannula.

Note: if a reusable obturator is used, remove it together with the guide.

Removal

- Insert the obturator through the cannula.
- While exerting continuous force upwards, turn the assembly counterclockwise until it is completely removed from the incision.

10. MATERIAL SPECIFICATIONS

Fergus cannula models FG-CN910670, FG-CN910690, FG-CN910870, and FG-CN910890 are made of the following materials:

- Obturator: acrylonitrile-butadiene-styrene (ABS)
- Seals: silicone
- Lateral portal: polyethylene (PE)
- Cannula: polycarbonate (PC)


Fergus cannula models FG-CN950670, FG-CN950690, FG-CN950870, and FG-CN950890 are made of the following materials:


- Obturator: acrylonitrile-butadiene-styrene (ABS)
- Seals: silicone
- Lateral portal: polyoxymethylene (POM)
- Cannula: polycarbonate (PC)

11. LEARN MORE


To learn more about this product, visit our website: www.fergusmedical.com.

12. SYMBOLS USED IN LABELS


	CATALOG NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	PRECAUTIONS

EXPIRATION DATE


DO NOT REUSE


MANUFACTURING DATE

MANUFACTURER


CHECK INSTRUCTIONS FOR USE

DO NOT USE IF PACKAGING IS DAMAGED


DO NOT RE-STERILIZE

KEEP AWAY FROM SUNLIGHT

KEEP DRY

**STERILE PRIMARY PACKAGING. STERILIZATION
METHOD: ETHYLENE OXIDE**

Fergus

PROMEDON S.A

Av. Gral. Manuel Savio s/n

Lote 3 • Manzana 3 • (X5123XAD)

Pque Ind. Ferreyra, Córdoba, Argentina

www.fergusmedical.com