Instrumentales no estériles

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Instructions for use
ENGLISH

M-350-00623 (Rev. 01)/11.07.22

Fergus

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1. DESCRIPTION OF THE DEVICE

The devices that conform the Fergus non-sterile instrument system may be single-use or reusable instruments.

These may be instruments used for resectioning, handling, piercing and cutting tissue, retrieving loose objects and handling sutures. The instruments have an ergonomic, low-profile design intended to gain manual control from different positions and for different functions.

The system includes instruments intended to handle sutures, such as suture passing instruments, which may be used to hold, stitch and retrieve sutures in a single effective step during an arthroscopic procedure. Other products used for cutting tissue and handling and managing suture are the suture-cutting and hinged instruments.

The instruments are designed to simplify surgical procedures for shoulder, knee, foot and ankle repair.

2. INDICATIONS FOR USE

The Fergus non-sterile instruments are indicated for shoulder, knee, foot and ankle arthroscopic surgical procedures, minimally invasive orthopedic surgery and/or open orthopedic surgery.

3. CONTRAINDICATIONS

- No contraindications were found for these devices.

4. WARNINGS

- Read the instructions for use carefully and completely before use.
- The device must be used only by a medical professional trained in the surgical techniques indicated.
- The law allows this device to be sold exclusively to physicians or people authorized by a physician.
- If the patient suffers or is suspected to suffer from bovine spongiform encephalopathy (BSE), products must be discarded as per local government and hospital procedures.
- The medical center is responsible for cleaning, sanitizing and sterilizing medical devices before use by following the local government regulatory requirements and these instructions for use.
- The instruments that have been designed for single use **MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED**, since this could impair their efficiency and increase the risk of inappropriate sterilization and cross-contamination.

- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative and/or governmental requirements.
- The device must not be used for surgical procedures that are not included in the "Indications for Use" section.
- These devices must not be implanted, since they are not intended for such use.
- Make sure that the product is disposed as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- Bio-hazardous waste, such as explanted devices, sutures, needles and contaminated surgical materials, must be disposed of safely and as per site policy.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the manufacturer.

5. PRECAUTIONS

- When two instruments are designed to be coupled, make sure that they are fully assembled before use.
- Excessive gripping or harsh manipulation may damage the locking mechanism of instruments with adjustable components.
- Do not use implant-specific instruments with other implants, since their effectiveness and safety cannot be guaranteed.
- The device must be removed before flexing the joint, since the device could bend or break.
- Handle sharp instruments with care to prevent possible injury.
- After inserting the instrument in the joint, excessive force or inappropriate movement could cause it to break and become lodged in soft tissue and/or disappear from the arthroscopic view of the surgical field, which may lead to fragments remaining inside the patient.
- Impacting instruments not suited for this process could cause damage to the device.
- Before using the device, examine it for damage. Do not use if damaged.
- When indicated, dismount the instrument to its component parts before cleaning, sanitization and sterilization.
- For strippers only: They must only be used to separate and remove soft tissue. They must not be used to protect soft tissue from saws or other surgical instruments.
- For needles for suture-passing forceps only: do not use outside the internal joint area, since it could cause damage to the joint.

- For suture-passing forceps with tip angled 20° upward and suture retrieval forceps with penetrating tip only:
 - The tip of the device must not be used as a lever or pivot point against bone or other hard tissues.
 - If the tip gets stuck, it must be removed by pulling the instrument backwards in a straight line. Do not turn, rotate or move the tip back and forth, since this could cause breakage.
 - Keep the jaws closed during penetration and open them only when you are ready to hold the desired suture.
 - Do not use the device to open tissue or apply excessive pressure on the device, since it could become damaged.
- For suture-cutting instruments only:
 - The safety mechanism must be unlocked before cutting the suture.
 - When the safety mechanism is locked, pressing the handle excessively could cause it to break rendering the instrument useless.
 - Do not use the device to open tissue or apply excessive pressure on the device, since it could become damaged.
 - Do not cut the tied suture knot if you are not seeing it directly, since it could become loose.
- For suture retrieval forceps only:
 - Do not use the device to penetrate, open or manipulate tissue, or apply excessive pressure on the device, since it could become damaged.
- For 4-mm spade-tip pins only:
 - Before passing the device, check that the knee flexion angles are appropriate.
 - The leg holder needs to be placed in the higher area of the leg, so that the pin outlet is exposed enough.
 - Do not bend the pin for insertion or during it, since this could damage the instrument.
- For sharp instruments, keep away from neurovascular structures.
- Medical centers must validate the cleaning, sanitization and sterilization cycles so as to guarantee the safety and the effectiveness of the instruments.

6. CLEANING AND SANITIZATION

Instruments must be cleaned, sanitized and sterilized before each use or reuse.

Single-use instruments must never be reused. A product is considered used after being opened, exposed to a non-sterile environment and/or coming into contact with blood, bone, tissue or other body fluids.

Reuse does not guarantee instrument safety and effectiveness, and it could create risks for user and/or patient health and/or safety.

Additionally, reusable instruments must be reprocessed, at most, two hours after use, by following explicitly the manufacturer indications for use and for preparing the cleaning solution.

Finally, products must be stored and transported inside a closed container to guarantee user safety.

The performance of these instruments is not affected when appropriately cleaned, sanitized and sterilized as mentioned in this document.

6.1. Detergent Selection

To select an appropriate detergent, consider the following:

- Detergent must be suitable for cleaning plastic and/or metallic biomedical instruments.
- It must be suitable for cleaning in ultrasonic machines (no foam formation).

The manufacturer recommends using enzyme or neutral-pH cleaners. Alkaline cleaning agents may be used if required by local laws.

Follow carefully the instructions of use provided by the detergent manufacturer and pay attention to the information on temperature of use, concentration, neutralization and rinsing.

6.2. Preliminary Cleaning

1. Take apart those instruments that can be disassembled before starting with the cleaning process.
2. Remove the dry visible residue from the instruments using cold water (under 40 °C). A sponge or a soft brush may be used. Pay special attention to joints, cavities, slits and other difficult to reach areas.
3. Rinse cavities, slits, joints and other difficult to reach areas and spaces energetically and carefully for over 30 seconds using a pressure water gun or clean tap water (with similar pressure).
4. Check that the devices have no visible traces of residue. If more residue is observed, repeat the previous steps until no more can be seen. If there is no visible residue on the products, proceed with manual cleaning and sanitization, or with machine (automatic) cleaning and thermal sanitization.

6.3. Manual Cleaning and Sanitization

1. Submerge the instruments in a solution of enzyme or alkaline detergent making sure that the instruments are opened and uncoupled.

Solutions with high alkalinity or low acidity are not recommended, since metallic or plastic parts could be damaged. If the solution does not have a neutral pH, rinsing, by means of validation at the end user site, and neutralization

must be guaranteed in order to ensure that the performance, fit and finish of the devices is appropriate.

The cleaning concentration used must follow detergent manufacturer specifications.

If this step can be done using an ultrasonic bath, do so for 10 minutes at a frequency of 35 kHz.

2. Rub the instruments for over 30 seconds using a soft brush to avoid damaging the surfaces, and pay special attention to those areas that may contain accumulated residue due to their geometry.

3. Rinse the devices carefully with distilled water for, at least, one minute.

4. Submerge the devices in a sanitizing solution, prepared following manufacturer instructions, for as long as the solution manufacturer establishes. Solutions with high alkalinity or low acidity are not recommended, since metallic or plastic parts could be damaged. If the solution does not have a neutral pH, rinsing, by means of validation at the end user site, and neutralization must be guaranteed in order to ensure that the performance, fit and finish of the devices is appropriate.

For hinged devices, operate the moving parts, at least, 5 times during the sanitization.

5. Rinse the instruments with distilled water or deionized sterile water for, at least, 1 minute.

6. Dry the instruments using medical grade compressed air.

7. Check that the devices have no visible traces of residue. If more residue is observed, repeat the previous steps until no more can be seen.

7. INSPECTION AND MAINTENANCE

- Fergus non-sterile instruments are precision devices and must be handled with care.
- The instruments must be inspected for damage before use.
- Cutting or sharp instruments could become dull with continuous use. This is not a defect, but rather normal device wear. Initial inspection must include the verification of the sharpness and cutting capability of these devices.
- Functional test of the manual operation of the instruments with an opening and closing mechanism: Check, without using excessive force, that the jaw and the tip are straight with no visible space when closed. Separate the handles and check that the movement is free and unrestricted, and that there is movement on the distal end. Repeat this process, at least, twice.
- Visual integrity inspection of the distal end: check

that the metal shaft of the flexible instruments is not cracked or fractured, since this indicates the end of life of the device.

- If any instrument presents damage, it must not be used. Instructions must be requested from the manufacturer for disposal and/or return.
- All mobile parts must be lubricated using biocompatible oil suitable for medical-grade instruments. Apply the lubricant following manufacturer recommendations.
- Perform a visual cleanliness inspection, checking for remaining residue.

8. PACKAGING

The devices must be packaged by using bags or wrapping material that meet the recommended specifications for steam sterilization in an autoclave. The packaging must be according to the ISO 11607 and EN 868 standards.

9. STERILIZATION

Before using any instruments that are not supplied sterile, you must clean them, sanitize them, place them in their final packaging and, finally, sterilize them.

Instruments must be steam-sterilized in an autoclave by using a fractional pre-vacuum process as per the ISO 17665-1 standard and considering the specific national requirements.

Sterilization must be performed by means of a fractional pre-vacuum process with the following parameters:

- 121 °C, time in vacuum of 20 minutes and 2 pre-vacuum cycles.

The operation instructions supplied by the autoclave manufacturer and the recommended guidelines for the maximum load of products to be sterilized must be followed, since they can change from one sterilizer to another. The autoclave must be validated, calibrated and maintained as per current local regulations.

Drying times may vary depending on the load used, and must increase as the load to be sterilized increases.

After sterilization, instruments must cool down without being touched or placed in a cold surface or submerged in cold liquid.

10. STORAGE

Non-sterile devices must be stored in a clean, dry and debris-free environment, protected from damage and at room temperature. The lifetime of these instruments is contingent to them being used appropriately as stated in these instructions.

The final user is responsible for the devices being stored after sterilization in a way that sterility is preserved until their use. If the sterility packaging is damaged, altered or shows signs of exposure to humidity, the instruments must be cleaned, packaged and sterilized again.

11. MORE INFORMATION

To learn more about this product, visit our website: www.fergusmedical.com

12. SYMBOLS USED IN LABELS

CATALOG NUMBER

BATCH NUMBER

PRECAUTIONS

MANUFACTURING DATE

NON-STERILE

MANUFACTURER

CHECK INSTRUCTIONS FOR USE

KEEP AWAY FROM SUNLIGHT

KEEP DRY

Fergus

PROMEDON SA
Av. Gral. Salvo s/n
Lote 3 • Manzana 3 • (X5123XAD)
Pque Ind. Ferreyra, Córdoba, Argentina

www.fergusmedical.com